

PRODUCT GROUP 9

HAND BRAKES, BRAKE HANDLES & RELEASE HANDLES




CONTENT

HAND BRAKES	204
Hand brakes - aluminium.....	204
Hand brakes - steel	205
Extension levers.....	206
Mounting clamps.....	207
BRAKE AND RELEASE HANDLES	208
Brake handles	208
Brake and release handles (double function).....	210
Release handles	211

HAND BRAKES

Hand brakes - aluminium


B1 - Brake with push to brake function and single bend brake lever



Description	Part. no. right	Part. no. left
Black	673061-R1	673062-L1

» Spring-loaded


B2 - Brake with push to brake function and double bend brake lever



Description	Part. no. right	Part. no. left
Black	673063-R1	673064-L1

» Spring-loaded

B3 - Brake with push to brake function and single bend brake lever



Description	Part. no. right	Part. no. left
Black	673065-R1	673066-L1

» Spring-loaded

B4 - Brake with push to brake function and double bend brake lever



Description	Part. no. right	Part. no. left
Black	673067-R1	673068-L1

» Spring-loaded

B5 - Brake with pull to brake function and single bend brake lever




Description	Part. no. right	Part. no. left
Black	673069-R1	673070-L1

» Spring-loaded

HAND BRAKES


B6 - Brake with pull to brake function and double bend brake lever



Description	Part. no. right	Part. no. left
Black	673071-R1	673072-L1

» Spring-loaded


B7 - Brake for sport with pull to brake function



Description	Part. no. right	Part. no. left
Black	673073-R1	673074-L1

» Spring-loaded

B8 - Brake for people with limited strength with push to brake function and double bend brake lever




Description	Part. no. right	Part. no. left
Black	673075-R1	673076-L1

» Spring-loaded

Hand brakes - steel

B9 - Brake with push to brake function and single bend brake lever



Description	Part. no. right	Part. no. left
Black	673077-R1	673078-L1

» Spring-loaded

HAND BRAKES

B10 - Brake with push to brake function and double bend brake lever



Description	Part. no. right	Part. no. left
Black	673079-R1	673080-L1


» Spring-loaded

Extension levers



Extension levers for MBL hand brakes

- » Ergonomic shape
- » Foldable which gives easy transfer in and out of the chair
- » Reduce the needed brake force on the brake handle
- » Available in 2 different lengths
- » Available for brakes with both 4mm and 6mm lever.
- » Nice elegant design
- » Material: Black PA+GF



Description	Part. no. right side	Part. no. left side
Extension lever for alu brakes (Short)	673219	673216
Extension lever for alu brakes (Long)	67321A	673217
Extension lever for steel brakes (Short)	67321G	673213
Extension lever for steel brakes (Long)	673218	673214



< Brake handle with short extension lever

Brake handle with long extension lever >




HAND BRAKES

Mounting clamps

Mounting clamp - black aluminium



Description	Part. no. Black
Complete mounting clamp for fitting brake on tube - Ø19/Ø20 mm	674005
Complete mounting clamp for fitting brake on tube - Ø19/Ø22 mm	674004
Complete mounting clamp for fitting brake on tube - Ø19/Ø25 mm	674001
Complete mounting clamp for fitting brake on tube - Ø19/Ø30 mm	674003



CUSTOMER SPECIFIED HAND BRAKES

Please note that the hand brakes shown represent our standard programme.

Should you have any interest in having hand brakes made after your specifications please do not hesitate to contact our sales personnel.


BRAKE AND RELEASE HANDLES

Automatic parking brake function




The stop lever pushes down automatically

Manual parking brake function



The stop lever pushes up - it has to be pushed down manually

Straight cable outlet vs. angled cable outlet



Angled cable outlet

Straight cable outlet

Brake handles



BH1 - Brake handle with angled cable outlet

- » Material: PA+GF
- » Textured surface

For tube diameter mm	Parking brake function	Quantity of parking brake positions	Type of screw for fitting	Color handle/stop lever	Part no.
Ø22	Automatic	5	Self tapping	Black/Black	671022
Ø22	Manual	5	Self tapping	Black/Black	671032
Ø25	Automatic	5	Self tapping	Black/Black	671009
Ø25	Manual	5	Self tapping	Black/Black	671019

BRAKE AND RELEASE HANDLES



BH1 - Brake handle with straight cable outlet

- » Material: PA+GF
- » Textured surface

For tube diameter mm	Parking brake function	Quantity of parking brake positions	Type of screw for fitting	Color handle/stop lever	Part no.
Ø22	Automatic	5	Self tapping	Black/Black	671020
Ø22	Manual	5	Self tapping	Black/Black	671030
Ø25	Automatic	5	Self tapping	Black/Black	671021
Ø25	Manual	5	Self tapping	Black/Black	671031



BH2 - Brake handles with angled cable outlet

- » Material: PA+GF
- » Textured surface
- » Version with automatic brake function for Ø25mm tube is available on request

For tube diameter mm	Parking brake function	Quantity of parking brake positions	Type of screw for fitting	Color handle/stop lever	Part no.
Ø22	Manual	5	Normal	Black/Red	671026
Ø22	Automatic	5	Self tapping	Black/Red	671029
Ø25	Manual	5	Normal	Black/Red	671027



BH2 - Brake handles with straight cable outlet

- » Material: PA+GF
- » Textured surface
- » Version with automatic brake function for Ø25mm tube is available on request

For tube diameter mm	Parking brake function	Quantity of parking brake positions	Type of screw for fitting	Color handle/stop lever	Part no.
Ø22	Manual	5	Self tapping	Black/Red	671042
Ø22	Automatic	5	Self tapping	Black/Red	671039
Ø25	Manual	5	Self tapping	Black/Red	671043

BRAKE AND RELEASE HANDLES

Brake and release handles (double function)



BHRH1 - Brake and release handle (double function) with angled cable outlet

- » Material: PA+GF
- » Textured surface
- » Available for Ø25 mm tube on request

For tube diameter mm	Parking brake function	Quantity of parking brake positions	Type of screw for fitting	Color handle/levers	Part no.
Ø22	Manual	5	Self tapping	Black/Black	671028



BHRH1 - Brake and release handle (double function) with straight cable outlet

- » Material: PA+GF
- » Textured surface
- » Available for Ø25 mm tube on request

For tube diameter mm	Parking brake function	Quantity of parking brake positions	Type of screw for fitting	Color handle/levers	Part no.
Ø22	Manual	5	Self tapping	Black/Black	671041

BRAKE AND RELEASE HANDLES

Release handles



RH1 - Release handle

- » Material: PA+GF
- » Textured surface
- » Available for Ø25 mm tube on request

For tube diameter mm	Type of screw for fitting	Color handle/lever	Part no.
Ø22	Self tapping	Black/Black	671034



RH2 - Release handle

- » Material: PA+GF
- » Textured surface
- » Available for Ø25 mm tube on request

For tube diameter mm	Type of screw for fitting	Color handle/lever	Part no.
Ø22	Normal	Black/Red	671044



RH3 - Release handle

- » Material: PA+GF
- » Textured surface
- » Available for Ø25 mm tube on request

For tube diameter mm	Type of screw for fitting	Color handle/lever	Part no.
Ø22	Normal	Black/Blue	671045

BRAKE AND RELEASE HANDLES



RH4 - Release handle

- » Material: PA+GF
- » Textured surface
- » Available for Ø25 mm tube on request

For tube diameter mm	Type of screw for fitting	Color handle/lever	Part no.
Ø22	Normal	Black/Black	67101E



RH5 - Release handle

- » Material handle: PA+GF
- » Material "holder": aluminium
- » Textured surface
- » Recommended for tube diameter Ø22-Ø27mm

For tube diameter mm	Color handle/holder	Part no.
Ø22	Black/Black	671105



Product group 9: Hand brakes, brake handles & release handles