


Handrims Titanium


Titanium – brushed – 41 mm long splices



- Very durable, yet lightweight
- Brushed, luxury looking surface
- 6 splices (D39005), 2 mm thick
- Splices can be adjusted +/- 2 mm
- Ø19 mm

Wheel Size	SD fitting distance mm	Weight g	Fitting	Part. no. Brushed
22"	466	370	6 splices	352282
24"	510	404	6 splices	352493
24"	517	409	6 splices	352482
25"	535	423	6 splices	352582
26"	567	444	6 splices	352682

Titanium – brushed – 37 mm long splices



- Very durable, yet lightweight
- Brushed, luxury looking surface
- 6 splices (D39016), 2 mm thick
- Splices can be adjusted +/- 2 mm
- Ø19 mm

Wheel Size	SD fitting distance mm	Weight g	Fitting	Part. no. Brushed
24"	510	400	6 splices	352449 (NEW)