

MBL line Front wheels 4"-8"




mb.l

Front wheels

4"–8"

Wheels with 3 spokes and solid tyre	4
Wheels with 5 spokes and solid tyre	7
Wheels with 7 spokes and solid tyre	11
Wheels with closed rim and solid tyre	13
Wheels with rubber tyre and tube or Flexel™ Lite S	16
Spacer bushes	17



All of our wheels are tested to the highest standards in our own testing laboratories.

Wheels with 3 spokes and solid tyre

4" - 8"

4" (100 x 26 mm) wheel with narrow, round, solid rubber or PU tyre



mm	mm	kg	RUBBER		PU	
			Grey	Black	Grey	Black
ø10	26	40	141406		15040D	15040E
ø8	36	40	140406		15040L	On request
ø8	45	40	140408		On request	15040I

4" (100 x 32 mm) wheel with wide, round, solid rubber or PU tyre



mm	mm	kg	RUBBER		PU	
			Grey	Black	Grey	Black
ø10	32	40	141426	141426B	15040G	15040H
ø8	34	40	140425	140425B	On request	150425B
ø8	36	40	140426	140426B	On request	On request
ø8	39	40	140498	On request	150498	150498S
ø8	40	40	140427	140427B	On request	On request
ø8	45	40	140428	140428B	on request	150428B




5" (125 x 27 mm) wheel with narrow, round, solid rubber or PU tyre



mm	mm	kg	RUBBER		PU	
			Grey	Black	Grey	Black
ø10	32	40	141504	141504B	On request	On request
ø8	36	40	140576	140576B	150576	On request
ø8	39	40	140598	On request	150598	150598S
ø8	40	40	140577	140577B	On request	On request




5" (125 x 34 mm) wheel with flat solid softroller rubber or PU tyre



 mm	 mm	 kg	RUBBER Grey	PU Grey
ø10	32	50	141572	on request
ø8	36	50	140571	150517




5.5" (140 x 40 mm) wheel with round solid softroller rubber or PU tyre



 mm	 mm	 kg	RUBBER Grey
ø10	41	50	141632
ø8	45	50	140632
ø8	53	50	140636
ø8	60	50	140634






6" (150 x 27 mm) wheel with narrow, round solid rubber or PU tyre



 mm	 mm	 kg	RUBBER Grey	PU Grey
ø10	32	40	141636	on request
ø8	45	40	140637	150637

6" (150 x 30 mm) wheel with ribbed polyurethane tyre






 mm	 mm	 mm	 mm	 kg	PU Grey
150	30	ø10	32	40	150652
150	30	ø8	36	40	150611
150	30	ø8	45	40	150613

Other dimensions are available on request.




6" (150 x 33 mm) wheel with wide, flat solid rubber or PU tyre



 mm	 mm	 kg	RUBBER Grey	PU Grey
ø10	38	40	141672	On request
ø8	45	40	140673	150673
ø8	60	40	140675	on request



7" (175 x 45 mm) wheel with ribbed polyurethane tyre



 mm	 mm	 kg	PU Grey
ø8	53	50	150732
ø8	60	50	150736

8" (200 x 50 mm) wheel with ribbed polyurethane tyre



 mm	 mm	 kg	PU Grey
ø8	53	50	150832
ø8	60	50	150836

Wheels with 5 spokes and solid tyre 5"-8"

5" (125 x 27 mm) wheel with round solid rubber or PU tyre



mm	mm	kg	RUBBER		PU	
			Grey	Black	Grey	Black
ø10	32	40	141506	141506B	150500	150506B
ø8	45	40	140508	140508B	On request	On request
ø8	53	40	140556	140556B	On request	On request

5" (125 x 40 mm) wheel with wide, flat solid rubber or PU tyre



mm	mm	kg	RUBBER		PU	
			Grey	Black	Grey	Black
ø10	41	40	141623	141623B	15061U	15061W
ø8	45	40	140623	140623B	15061C	15061B
ø8	60	40	140625	140625B	150525	150525B




5.5" (140 x 40 mm) wheel with black round solid softroller rubber or PU tyre



mm	mm	kg	RUBBER		PU	
			Grey	Black	Grey	Black
ø10	10	50	141622	141622B	15061K	15061I
ø8	45	50	140622	140622B	150618	150618B
ø8	60	50	140624	140624B	15060Z	15061D

6" (150 x 30 mm) wheel with ribbed polyurethane tyre



 mm	 mm	 kg	PU Grey
ø10	32	40	150651

7" (175 x 43 mm) wheel with smooth polyurethane tyre



 mm	 mm	 kg	PU Grey
ø8	53	50	150719
ø8	60	50	150720




Sign up to our
newsletter at
mbl.dk to stay
updated!








7" (175 x 45 mm) wheel with ribbed polyurethane tyre



 mm	 mm	 kg	PU Grey
ø8	60	50	150726


8" (200 x 45 mm) wheel with smooth polyurethane tyre



 mm	 mm	 kg	PU Grey
ø8	60	50	150820

8" (200 x 50 mm) wheel with ribbed polyurethane tyre






 mm	 mm	 kg	PU Grey
ø8	53	50	150822
ø6	56	50	150824
ø8	60	50	150826

Wheels with 7 spokes and solid tyre 6"-8"




6" (150 x 27 mm) wheel with round solid rubber or PU tyre




 mm	 mm	 kg	RUBBER Grey	PU Grey
ø10	32	75	141606	on request
ø8	36	75	140606	on request
ø8	40	75	140655	150655
ø8	45	75	140608	150610




6" (150 x 33 mm) wheel with flat solid rubber or PU tyre



 mm	 mm	 kg	RUBBER Grey	PU Grey
ø10	32	75	141616	15061J
ø8	36	75	on request	150620
ø8	40	75	on request	on request
ø8	45	75	140618	on request
ø8	53	75	on request	on request
ø8	56	75	on request	on request
ø8	58	75	-	on request
ø8	60	75	140668	on request




7" (175 x 27 mm) wheel with round solid rubber or PU tyre



 mm	 mm	 kg	RUBBER Grey	PU Grey
ø10	37	40	141706	on request
ø8	45	40	140706	150708
ø8	53	40	140756	on request

7" (175 x 33 mm) wheel with flat solid rubber or PU tyre



 mm	 mm	 kg	RUBBER Grey	PU Grey
ø10	37	50	141716	150707
ø8	45	50	140716	150705




7½" (190 x 29 mm) wheel with narrow, smooth/ribbed polyurethane tyre, different color combinations




 mm	 mm	 kg	Color rim	PU Black	PU Grey
ø8	53	40	Grey	on request	150744




8" (200 x 27 mm) wheel with round solid rubber or PU tyre



 mm	 mm	 kg	RUBBER Grey	PU Grey
ø10	37	40	141806	15080Z
ø8	45	40	140806	150810




8" (200 x 33 mm) wheel with flat solid rubber or PU tyre



 mm	 mm	 kg	RUBBER Grey	PU Grey
ø10	37	40	141816	on request
ø8	60	40	140818	15081F




8" (200 x 40 mm) wheel with 7 spokes and smooth V-shaped polyurethane tyre



 mm	 mm	 kg	PU Grey
ø8	60	40	15080L

8" (200 x 30 mm) wheel with narrow, ribbed polyurethane tyre





 mm	 mm	 kg	PU Grey
ø8	60	40	150805

Wheels with closed rim and solid tyre 7"-8"




8" (200 x 50 mm) wheel with rim in one piece and ribbed polyurethane tyre



 mm	 mm	 kg	PU Grey
ø8	60	50	150855

8" (200 x 50 mm) wheel with split rim and ribbed polyurethane tyre



 mm	 mm	 kg	PU Grey
ø8	53	60	150812

Carbon solutions for rehab industry

**Skilfully engineered. Delivered fast.
Quality you can trust.**






View the video about
our carbon competencies.

Wheels with rubber tyre and tube or Flexel™ Lite S 8"






8" (200 x 50 mm) wheel with grey tyre




 mm	 mm	 kg	Pneumatic Grey	FlexelLite Grey
ø8	53	50	160832	on request
ø8	60	50	160836	on request




8" (200 x 50 mm) wheel with grey tyre




 mm	 mm	 kg	Pneumatic Grey	FlexelLite Grey
ø8	60	50	160826	on request

8" (200 x 30 mm) wheel with grey tyre



 mm	 mm	 kg	Pneumatic Grey	FlexelLite Grey
ø10	37	40	160852	on request
ø6	56	40	160807	on request

Spacer bushes



- » To match the correct hub width to fork width
- » For ø10 mm axles
- » Part. numbers are for packings with 200 pcs
- » Minimum order is 200 pcs. per order



Choice of spacer bushes for all hub widths:

mm	26 mm			32 mm			37 mm		
	B	D = 6	D = 8	B	D = 6	D = 8	B	D = 6	D = 8
36	5	19406E	19408E	2	19406B	19408B	-	-	-
37							-	-	-
38	6	19406F	19408F	3	On request	19408C		-	-
39							1	19406A	19408A
40	7	19406G	19408G	4	19406D	19408D			
41							2	19406B	19408B
42	8	On request	19408H	5	On request	19408E			
43							3	19406C	19408C
44	9	19406I	19408I	6	On request	19408F			
45							4	19406D	19408D
46	10	On request	19408J	7	19406G	19408G			
47							5	19406E	19408E
48	11			8					
49							6	19406F	19408F
50	12		19408L	9	19406I	19408I			
51							7	19406G	19408G
52	13	On request	-	10	19406J	19408J			
53							8	19406H	19408H
54	14	On request	19408N	11	On request	19408K			
55		-	-				9		19408I
56		-	-	12	On request	19408L	9,5	*19406I	
67		-	-				10	19406J	19408J
58		-	-	13	On request	-			
59		-	-				11	19406K	19408K
60		-	-	14	On request	19408N			
61		-	-				12	19406L	19408L

Would you like to see and learn more?

MBL website



Our videos on YouTube



Check our Facebook



Connect with us on LinkedIn



Connect with us on WeChat



Other product categories:

Rear wheels 12"-26"



Tyres, PU tyres, tubes , infills & tyres
w insert for manual wheelchairs 6"-26"



Pushrims & pushrim covers



Front forks & bearing housings



QR axles and bushings



Hand brakes, brake handles & release handles



Other wheelchair components & accessories



Foam parts, supports & textiles



Tyres, tubes and tyres with infills
for power wheelchairs



MBL Denmark

MBL a/s

Glarmestervej 18 B, St. Th.

8600 Silkeborg

CVR: 12825242

contact@mbl.dk

MBL Poland

MBL Poland Sp. z o.o.

Sulejowska 45d

97-300 Piotrkow Trybunalski

VAT No: 7710008003

Tel: +48 44 73 24 100

contact@mbl.pl

MBL China

MBL (Xiamen) Co. Ltd.

No.289 Wengjiao Road Xinyang

Industrial District Haicang

361022 Xiamen

Tel: +86 592 5961 896

contact@mbl.com.cn

March 2026