

Tyres for manual wheelchairs

Air pressure explained

1 Psi = 0,069 bar Bar is the atmospheric pressure at the sea level.
 1 bar = 14,5 Psi Psi is 1 pound of force per square inch.


ETRTO explained

An ETRTO size of 47-507 means that the tyre is 47 mm wide and has an inner diameter of 507 mm.

Tyres for rear wheels 20"-26"




Size	ETRTO	Pressure psi/bar	Part no. grey	Part no. black
24" x 1.75" x 2"	47-507	90/6	on request	432436




Size	ETRTO	Pressure psi/bar	Part no. grey	Part no. black
20" x 1 3/8"	37-451	65/4.5	412061	on request
22" x 1 3/8"	37-489	65/4.5	412261	432261
22" x 1 3/8"	37-501	65/4.5	412263	on request
24" x 1 3/8"	37-540	65/4.5	412461	41240F

NON MARKING



Size	ETRTO	Pressure psi/bar	Part no. grey	Part no. black
22" x 1 3/8"	37-501	65/4.5	412278	on request
24" x 1 3/8"	37-540	65/4.5	412478	432401N





Size	ETRTO	Pressure psi/bar	Part no. grey	Part no. grey/black
20" x 1"	25-451	110/7.5	412051	412054
22" x 1"	25-489	110/7.5	412251	
22" x 1"	25-498	110/7.5	412252	412255
24" x 1"	25-540	110/7.5	412451	412454
25" x 1"	25-559	110/7.5	412551	
26" x 1"	25-590	110/7.5	412651	



Size	ETRTO	Pressure psi/bar	Weight approx. g	Part no. grey	Part no. grey/black
24" x 1"	25-540	110/7.5	410	412471	412474


NON MARKING is an option for all our tyres.

MOQ 300-500 pieces


Please contact our sales department.



Tyres for wheels 6" - 12 1/2"




Size	ETRTO	Pressure psi/bar	Part no. grey	Part no. black
12 1/2" x 2 1/4"	57-203	35/2.5	411262	431262
16" x 1.75"	47-305	35/2.5	411662	431662



Size	ETRTO	Pressure psi/bar	Part no. grey	Part no. black
7" x 1 3/4"	47-94	36/2.5	410751	on request



Size	ETRTO	Pressure psi/bar	Part no. grey	Part no. black
6" x 1¼"	32-87	36/2.5	410652	-
8" x 1¼"	32-137	36/2.5	410852	on request
200 x 50	50-94	36/2.5	410853	on request



Size	ETRTO	Pressure psi/bar	Part no. grey	Part no. black
12½" x 2¼"	57-203	30/2.1	411271	on request




Size	ETRTO	Pressure psi/bar	Part no. grey	Part no. black
12½" x 2¼"	57-203	30/2.1	411254	on request

Sign up to our
newsletter at
mbl.dk to stay
updated!




About valves

Valves - two types:




Car valve/Schrader




Racer valve/Presta


Valves - varied shapes:




Alfa




Beta




α/β



90/90



45/45



90/30

View from the side

View from above

Tubes for manual wheelchairs

5"-26"

Size	ETRTO	Valve angle /	Part no. Schrader/Car valve
6" x 1¼"	32-87	45/45	450601
8" x 1¼"	32-137	45/45	450805
8" x 2"	50-94	90/30	450802
200 x 50	50-94	90/90	450801
12½" x 2¼"	62-203	Straight	451202
12½" x 2¼"	62-203	90/90	451203
12½" x 2¼"	62-203	45/45	451213
16" X 1,75" X 2,125"	47/57-305	Straight	451601
20" x 1 3/8"	37-451	Straight	452001
20" x 1.75" x 2.125"	47/57-406	Straight	452002
20" x 1"	25-451	Straight	452005
22" x 1 3/8"	37-489/501	Straight	452201
22" x 1"	25-489/501	Straight	452202
24" x 1"	25-540	Straight	452406
24" x 1 3/8"	37-540	Straight	452403
25" x 1"	25-559	Straight	45250
25" x 1" (extra light)	18/25-559	Straight	-
26" x 1"	25-590	Straight	452601