

NOVA handrims with anti-slip insert

The NOVA handrim is a good choice for users with full or slightly reduced hand strenght. The anti-slip insert made of environmentally friendly food grade silicone provides an excellent grip which makes the wheelchair more maneurable and easier to drive.

- High sides of hard anodized aluminium for easier braking
- The silicone insert is not too sticky which results in a reliable and comfortable grip as sand and other materials do not stick onto it.
- Good anti-slip performance even when wet.
- 6 splices (39210Q)
- Available in black hard anodized – shiny silver anodized on request
- Replaceable silicone insert



Wheel Size	SD fitting distance mm	Weight g	Silver anodized (Shiny) NOVA S	Black hard anodized NOVA H
24"/507	472	-	-	34228S (NEW)
24"	510	570	coming	34248G (NEW)
24"	515	590	-	34248S
25"	530	615	coming	34258G (NEW)
25"	535	630	-	34258S
26"	565	690	-	34268S

