

CYCLONE

24"

CYCLONE is a strong yet elegant looking wheel. The slim hub and only 12 spokes give the wheel a light look and fitted with the Omobic SpeedLite tyre it will guarantee the user a very easy run. The hard anodized pushrim is durable and generates minimum heat when braking.






Black rim
(DW2)



Aluminium pushrim
6 splices with two holes
Black or silver anodized



- Black anodized hub (R6) 48/70mm
- 12 black KTL steel spokes Ø3 mm



Omobic SpeedLite
25 mm



Max load
60 kg

Raw wheels

SIZE	SD diameter mm	Weight g	Part no. ø12 mm	Part no. ø1/2"
24"/Ø540	510	-	234AD3	234AD2

Wheels with pushrim and tyre (directional)

SIZE	Pushrim	g	Part no. ø12 mm LEFT	Part no. ø12 mm RIGHT	Part no. ø1/2" LEFT	Part no. ø1/2" RIGHT
24"	Black anodized	1550	264AD8-L	264AD8-R	264AD2-L	264AD2-R
24"	Silver anodized	1550	264ADC-L	264ADC-R	264AD9-L	264AD9-R

Quick-release axles in steel or titanium

Material	Weight (g) Ø12 mm / Ø1/2"	Lenght mm	Release button colour	Part. no. Ø12 mm	Part. no. Ø1/2"
Steel	97/108	108	Black	618254	618354
Titanium	66/73	108	Black	61825A	61835A