

## DAYLITE

24"-25"

The DAYLITE is a beautiful and elegant wheel for everyday use and an absolute favorite among many customers! The complete 24" wheel weighs only 1060 g. The radial lightweight hub was designed by MBL's founder and is now a worldwide classic. Fitted with the OMOBIC lightweight yet durable SpeedLite tyre the DAYLITE is a light, durable and stylish addition to most wheelchairs.





Black rim  
(DW4)



Extra light aluminium  
pushrim  $\varnothing 19 \times 0,8$  mm.  
6 splices, one hole  
Black or silver anodized



Black hub (R1)  
70/48mm  
24 extra light stainless  
steel spokes



Omobic SpeedLite  
25 mm



Max load  
60 kg

## Raw wheels

SIZE	SD diameter mm		Part no. $\varnothing 12$ mm	Part no. $\varnothing 1/2$ "
24"	517	525	234L11	234L61
25"	536	540	235L11	235L61

## Wheels with pushrim and tyre

SIZE	Pushrim	Weight g	Part no. $\varnothing 12$ mm	Part no. $\varnothing 1/2$ "
24"	Silver anodized	1060	264L11	264L61
24"	Black anodized	1060	264L12	264L62
25"	Silver anodized	1080	265L11	265L61
25"	Black anodized	1080	265L12	265L62

## Quick-release axles in steel or titanium

Material	Weight (g) $\varnothing 12$ mm / $\varnothing 1/2$ "	Lenght mm	Release button colour	Part. no. $\varnothing 12$ mm	Part. no. $\varnothing 1/2$ "
Steel	97/108	108	Black	6180L1	6181L1
Titanium	66/73	108	Black	6180L8	6181L8