

LOCK-A wheel lock in aluminium



- Luxury and very efficient wheel-lock
- Easy and smooth operation – can be used also by users with reduced hand strength
- Easy to install on any existing clamp systems
- Ergonomic hand grip
- Fit on most chairs and suitable for all kind of wheelchairs 90 degrees angle

Description	Weight g	Part no.
Activate brake with LEFT hand	130	67303Z
Activate brake with RIGHT hand	130	67303Y

LOCK-P wheel lock in plastic



- Easy to install on any existing clamp systems
- Can be installed on both sides of wheelchair, to the users preference
- Ergonomic hand grip
- Easy and smooth operation
- Fit on most chairs and suitable for all kind of wheelchairs
- 90 degrees angle

Description	Weight g	Part no.
Activate brake with LEFT hand	100	67303MA
Activate brake with RIGHT hand	100	67303NA

Clamp



Fitting Tube Dimension	Part no.
Ø19/Ø20	674005
Ø19/Ø22	674004
Ø19/Ø25	674001
Ø19/Ø30	674003