

MAGNILITE

24"

The MAGNILITE wheel is built with our DW2 rim profile. The magnesium rim weighs only 365 g compared to 476 g for the same rim in aluminium and combined with our Micro Hub, 18 extra light spokes and a super light tyre you get a beautiful wheel with a weight just around 850 g for the complete wheel. This wheel is extremely light and easy to handle - it simply makes life on the go way easier!





Magnesium
(DW2)
Black or white



Extra light aluminium
pushrim Ø19 x 0,8 mm.
6 splices, one hole
Black or silver anodized



Black micro hub (R5),
48/70mm, 18 extra light
stainless steel spokes

18 extra light stainless steel
spokes



24" Omobic Gazelle
25 mm



Max load
40 kg

Raw wheels

RIM COLOR	SIZE	SD diameter mm	Weight g	Part no. ø12mm	Part no. ø1/2"
Black	24"	510	-	234MG1	234MG2
White	24"	510	-	On request	On request

Wheels with pushrim and tyre

RIM COLOR	SIZE	Pushrim	Weight g	Part no. ø12 mm	Part no. ø1/2"
Black	24"	Silver anodized	850	264MG11	264MG61
Black	24"	Black anodized	850	264MG12	264MG62
White	24"	Black anodized	850	264MG89	264MG99

Quick-release axles in titanium

Material	Weight (g) Ø12 mm / Ø1/2"	Lenght mm	Release button colour	Part. no. Ø12 mm	Part. no. Ø1/2"
Titanium	63/70	103	Black aluminium	6180SF	6181SF