


NOVA pushrims with anti-slip insert

The NOVA pushrim is a good choice for users with full or slightly reduced hand strength. The anti-slip insert made of environmentally friendly food grade silicone provides an excellent grip which makes the wheelchair more maneuverable and easier to drive.

- High sides of hard anodized aluminium for easier braking
- The silicone insert is not too sticky which results in a reliable and comfortable grip as sand and other materials do not stick onto it.
- Good anti-slip performance even when wet.
- 6 splices (39210Q)
- Available in black hard anodized – shiny silver anodized on request
- Replaceable silicone insert



| Wheel Size | SD fitting distance mm | Fit Omobic rim types | Weight g | Black hard anodized NOVA H |
|------------|------------------------|----------------------|----------|----------------------------|
| 24" | 515 | DW2 + DW4 + DW7 | 590 | 34248S |
| 24" | 523 | SW1 | 610 | 34248T |
| 25" | 535 | DW2 + DW4 + DW7 | 630 | 34258S |
| 26" | 565 | DW2 + DW4 + DW7 | 690 | 34268S |

