

## CHALLENGER with off-set hub

24"

The OFF-ROAD CHALLENGER is an optimal choice for active users who often come across bumpy and challenging roads. The OFF-ROAD CHALLENGER is a strong and durable wheel equipped with a 10 mm off-set high flange hub. The 28 mm tall rim profile results in a dynamic and sporty look.





Black rim  
High profile



Aluminium pushrim  
6 splices with two holes  
Black anodized  
Ø19 mm



10 mm off-set high flange  
hub (F5) 53/64 mm  
36 stainless steel spokes.



- Innova tyre  
IA2003, 47 mm



Max load  
60 kg

## Raw wheels

SIZE	SD diameter mm	Weight g	Part no. ø12 mm	Part no. ø1/2"
24"/Ø507	472	970	234EZA	234EZA

## Wheels with pushrim and tyre

SIZE	Tyre DIRECTIONAL	Weight g	Part no. ø12 mm	Part no. ø1/2"
24"/Ø507	LEFT	2120	264EXG-L	264EZG-L
24"/Ø507	RIGHT	2120	264EXG-R	264EZG-R

## Quick-release axles in steel or titanium

Material	Weight (g) Ø12 mm / Ø1/2"	Length mm	Release button colour	Part. no. Ø12 mm	Part. no. Ø1/2"
Steel	97/108	108	Black	618254	618354
Titanium	66/73	108	Black	61825A	61835A