


## TWIN – for power assisted wheelchairs

24"

Our OFF-ROAD TWIN WHEEL will take you anywhere! It is designed for use with power assisted wheelchair systems. The two tyres make a large surface, especially when a little "flat". Each tyre is 6 cm wide, so with two wheels in one wheel we get a surface that is 12 cm wide. This makes the wheel suitable for beaches etc. The TWIN wheel is most often used without pushrim and cambering of the wheels is important for stabilizing the wheelchair at high speed





Black rim (SW1)  
Two rims welded together



Aluminium pushrim  
6 splices with two holes  
Black anodized  
Ø19 mm



Black off-set hub 7mm, with high flanges - ekstra long 90/79 mm  
36 stainless steel spokes, 2 mm



Innova tyre IA2003,  
47 mm



Max load  
60 kg

## Wheels with tyre and pushrim

SIZE	Pushrim	Tyre DIRECTIONAL	Weight g	Part no. ø12 mm	Part no. ø1/2"
24"/Ø507	Black anodized aluminium	LEFT	-	-	264TWIN98L
24"/Ø507	Black anodized aluminium	RIGHT	-	-	264TWIN98R

## Wheels with tyre – no pushrim

SIZE	Pushrim	Tyre DIRECTIONAL	Weight g	Part no. ø12 mm	Part no. ø1/2"
24"/Ø507	-	LEFT	3200	-	244TWIN98L
24"/Ø507	-	RIGHT	3200	-	244TWIN98R

## Quick-release axles in steel

Lenght mm	Release button colour	Part. no. Ø12 mm	Part. no. Ø1/2"
130	Black	-	618326