


Optimum ZF with anti-slip coating



- Optimal for users with reduced hand strength
- Available with 6 splices (39210Q) or 6 rivnuts
- 1,5 mm anti-slip coating, Lizard grip

Surface treatment	Wheel Size	Type of fitting	SD / Rivnut Fitting distance	Fit Rim Type	Weight g	Item no. Produced in Poland
Optimum ZF - Anti-slip Coating	24"	6 splices	515	DW2 + DW4 + DW7	560	3424F5-D
Optimum ZF - Anti-slip Coating	24"	6 splices	523	SW1	530	3424J1-D
Optimum ZF - Anti-slip Coating	25"	6 splices	535	DW2 + DW4 + DW7	550	3425F5-D
Optimum ZF - Anti-slip Coating	26"	6 splices	565	DW2 + DW4 + DW7	560	3426F5-D
Optimum ZF - Anti-slip Coating	24"	6 rivnuts	515	SW1	535	32244M-D
Optimum ZF - Anti-slip Coating	24"	6 rivnuts	520	DW5	550	32244N-D

Optimum Z



- The large and ergonomic shape gives the user a good and comfortable grip
- Available with 6 splices (39210Q)
- Black hard anodized aluminium - very scratch resistant

Surface treatment	Wheel Size	SD Fitting distance	Fit Rim Type	Weight	Item no.
Optimum ZH - black hard anodized	24	515	DW2 + DW4 + DW7	275	3424F5
Optimum ZH - black hard anodized	24	523	SW1	280	3424F4
Optimum ZH - black hard anodized	25	535	DW2 + DW4 + DW7	285	3425F5
Optimum ZH - black hard anodized	26	565	DW2 + DW4 + DW7	300	3426F5

