



Optimum F with anti-slip coating



- Optimal for users with reduced hand strength
- Available with 6 splices (39210Q) or 6 rivnuts
- 1,5 mm anti-slip coating, Lizard grip

Surface treatment	Wheel Size	Type of fitting	SD / Rivnut Fitting distance	Fit Rim Type	Weight g	Item no. Produced in Poland
Optimum F - Anti-slip coating	24	6 splices	515	DW2 + DW4 + DW7	520	34244S-D
Optimum F - Anti-slip coating	24	6 splices	523	SW1	530	3424DS-D
Optimum F - Anti-slip coating	25	6 splices	535	DW2 + DW4 + DW7	540	34254S-D
Optimum F - Anti-slip coating	26	6 splices	565	DW2 + DW4 + DW7	560	34264S-D
Optimum F - Anti-slip coating	24	6 rivnuts	515	SW1	510	32244R-D
Optimum F - Anti-slip coating	24	6 rivnuts	520	DW5	520	32244P-D

Optimum



- The large and ergonomic shape gives the user a good and comfortable grip
- Available with 6 splices (39210Q) or 6 rivnuts
- Black hard anodized aluminium - very scratch resistant

Surface treatment	Wheel Size	Type of fitting	SD / Rivnut Fitting distance	Fit Rim Type	Weight	Item no.
Optimum H - Black hard anodized	24	6 splices	515	DW2 + DW4 + DW7	275	34244S
Optimum H - Black hard anodized	24	6 splices	523	SW1	280	3424DS
Optimum H - Black hard anodized	25	6 splices	535	DW2 + DW4 + DW7	290	34254S
Optimum H - Black hard anodized	26	6 splices	565	DW2 + DW4 + DW7	300	34264S
Optimum H - Black hard anodized	24	6 rivnuts	515	SW1	270	32244R

