

SKYLITE

20"-25"

SKYLITE is your ultimate choice when you want to travel elegant and light. Our engineers have, as the first in the world, developed and tested micro hubs for wheelchair use. Combined with an ultra lightweight rim, where we have used a special CNC milling method between the spoke holes, we have developed a complete wheel weighing only 960 grams.





Black rim with silver cavities to reduce weight (DW7)



Extra light aluminium pushrim $\varnothing 19 \times 0,8$ mm.
6 splices, one hole
Black or silver anodized



Polished micro hub (R5)
70/48 mm
18 extra light stainless steel spokes



20", 22", 25": OMOBIC SpeedLite, 25 mm



Max load
40 kg



24. OMOBIC Gazelle
25 mm. Gazelle is optional for the 25" wheel.

Raw wheels

SIZE	SD diameter mm	Weight g	Part no. $\varnothing 12$ mm	Part no. $\varnothing 1/2$ "
20"	428	410	230US1	230US2
22"	466	430	232US1	232US2
24"	517	455	234US1	234US2
25"	536	490	235US1	235US2

Wheels with pushrim and tyre

SIZE	Pushrim	Tyre	Weight g	Part no. $\varnothing 12$ mm	Part no. $\varnothing 1/2$ "
20"	Silver anodized	OMOBIC Speedlite	925	260US1	260US2
20"	Black anodized	OMOBIC Speedlite	925	260US3	260US4
22"	Silver anodized	OMOBIC Speedlite	955	262US1	262US2
22"	Black anodized	OMOBIC Speedlite	955	262US3	262US4
24"	Silver anodized	OMOBIC Gazelle	970	264US11	264US12
24"	Black anodized	OMOBIC Gazelle	970	264US13	264US14
25"	Silver anodized	OMOBIC Speedlite*	990	265US1	265US2
25"	Black anodized	OMOBIC Speedlite*	990	265US3	265US4

* Gazelle ultralight tyre is now available in 25" if you wish the lightest option of tyre

Quick-release axles in steel or titanium

Material	Weight (g) $\varnothing 12$ mm / $\varnothing 1/2$ "	Lenght mm	Release button colour	Part. no. $\varnothing 12$ mm	Part. no. $\varnothing 1/2$ "
Steel	92/103	103	Mirror polished	6180SD	6181SD
Titanium	63/70	103	Mirror polished	6180SA	6181SA