

SMARTLITE

24"-25"

The SMARTLITE is an elegant wheel with a single wall rim for everyday use. The radial lightweight hub was designed by our founder and is now a worldwide classic. Now fitted with the Omobic lightweight yet durable SpeedLite tyre the SMARTLITE will be a light and stylish addition to most wheelchairs (Weight with OMOBIC SpeedLite tyre and pushrim: 24"/1250 g)





Black rim
(SW1)



Extra light aluminium
pushrim Ø19 x 0,8 mm.
6 splices, one hole
Black or silver anodized



Black hub (R1)
70/48mm
24 extra light stainless
steel spokes



Omobic SpeedLite
25 mm



Max load
60 kg

Raw wheels

SIZE	SD diameter mm	Weight g	Part no. Ø12mm	Part no. Ø1/2"
24"	523	670	234ES11	234ES61
25"	542	700	235ES11	235ES61

Wheels with pushrim and tyre

SIZE	Pushrim	Weight g	Part no. Ø12 mm	Part no. Ø1/2"
24"	Silver anodized	1240	264ES11	264ES61
24"	Black anodized	1240	264ES12	264ES62
25"	Silver anodized	1300	265ES11	265ES61
25"	Black anodized	1300	265ES12	265ES62

Quick-release axles in steel or titanium

Material	Weight (g) Ø12 mm / Ø1/2"	Length mm	Release button colour	Part. no. Ø12 mm	Part. no. Ø1/2"
Steel	97/108	108	Black	6180L1	6181L1
Titanium	66/73	108	Black	6180L8	6181L8