


Pushrims Titanium

Titanium – brushed



- Very durable, yet lightweight
- Brushed, luxury looking surface
- 6 splices (D39005), 2 mm thick
- Splices can be adjusted +/- 2 mm
- Ø19 mm

Wheel Size	SD fitting distance mm	Fit Omobic rim types	Weight g	Fitting	Part. no. Brushed
22"	466	DW4 & DW7	370	6 splices	352282
24"	510	DW2	404	6 splices	352493
24"	517	DW4 & DW7	409	6 splices	352482
25"	535	DW4 & DW7	423	6 splices	352582
26"	567	DW4 & DW7	444	6 splices	352682



Looking for something special?

Ask sales for customized pushrims!