


TRIP

Folding handle



Features

- The push handle is foldable to save packing and storage space
- The handle has an ergonomic shape which makes it comfortable for the wheelchair assistant.
- TPR stripes on the top and the bottom of the handle provide a better grip and relieve hand pain and soreness
- The metal tube can be customized to fit specific dimensions of wheelchair back canes



A Tube diameter mm	B mm	C mm	H mm	L mm	Part no.
ø22,8	without	without	without	45	67325R
ø14,5	10,1	25	2 x M6	45	67325G
ø14,5	8	-	1 x M5	25	67325X
ø14,7	10,1	25	2 x M6	45	67325B
ø16	without	without	without	45	67325C
ø20	12	22	2 x M6	45	67325J
ø15	10	15	2 x M5	45	67325N
ø21	without	without	without	25	67325S