

# Compact and lightweight air inflator

## Its advantage over traditional solutions are many:

- The inflator measures only 234 x 57 x 45 mm
- Weight: 430 g
- Easy use, it inflates the air automatically. All you need is to press the button
- No battery requirement. The air inflator is charged from a socket (type-C input). Charging time 2 hours
- Up to 10BAR inflating
- Due to the fact that the set contains adapters for different types of valves, the air inflator has a wide range of purposes
- The above features of our inflator give the user the possibility to instantly inflate deflated wheelchair wheels and regain mobility and independence

Charging cable, inflating attachment with valve type and ball converters and users manual are included with the product.



Description	Part no.
Air inflator A10, Black, Max Pressure: 150 PSI, with LED display	EBE01D-O