


Aluminium pushrims

The originals

Aluminium



- Ø19 x 1,25
- Splices (392115) with two holes
- Extra long splices used for pushrims which fit the Buddy wheels
- Splices can be adjusted +/- 2 mm

Wheel Size	SD /rivnut fitting distance (mm)	Fit rim types	Weight g	Fitting	Splice Type	Silver anodized	Black anodized
20"	422	DW2	-	6 splices	392115	on request	on request
20"	434	SW1	292	6 splices	392115	342051	342052
22"	460	DW2	-	6 splices	392115	on request	on request
22"	466	DW4 & DW7	311	6 splices	392115	342255	342256
22"	472	SW1	315	6 splices	392115	342251	342252
24"	480	Buddy SW Wide	-	6 splices	392144	-	3424E4
24"	510	DW2	338	6 splices	392115	344458	342458
24"	517	DW4 & DW7	343	6 splices	392115	342455	342456
24"	523	SW1	350	6 splices	392115	342451	342452
25"	530	DW2	351	6 splices	392115	342568	342567
25"	535	DW4 & DW7	355	6 splices	392115	342555	342556
25"	542	SW1	359	6 splices	392115	342561	342558
26"	560	DW2	-	6 splices	392115	on request	on request
26"	567	DW4 & DW7	374	6 splices	392115	342655	342656
26"	573	SW1	378	6 splices	392115	342651	342652
22"	469	DW5	298	6 rivnuts	398013	322284	322285
24"	515	SW1	327	6 rivnuts	398013	322452	322450
24"	520	DW5	330	6 rivnuts	398013	322495	32240BX
25"	539	DW5	341	6 rivnuts	398013	322595	322594