



## Optimum F with anti-slip coating



- Optimal for users with reduced hand strength
- Available with 6 splices (39210Q) or 6 rivnuts
- 1,5 mm anti-slip coating, Lizard grip

| Surface treatment             | Wheel Size | Type of fitting | SD / Rivnut<br>Fitting distance | Weight<br>g | Item no.<br>Produced in Poland |
|-------------------------------|------------|-----------------|---------------------------------|-------------|--------------------------------|
| Optimum F - Anti-slip coating | 24         | 6 splices       | 515                             | 520         | 34244S-D                       |
| Optimum F - Anti-slip coating | 25         | 6 splices       | 535                             | 540         | 34254S-D                       |
| Optimum F - Anti-slip coating | 26         | 6 splices       | 565                             | 560         | 34264S-D                       |
| Optimum F - Anti-slip coating | 24         | 6 rivnuts       | 515                             | 510         | 32244R-D                       |

## Optimum – black hard anodized



- The large and ergonomic shape gives the user a good and comfortable grip
- Available with 6 splices (39210Q) or 6 rivnuts
- Black hard anodized aluminium - very scratch resistant

| Surface treatment               | Wheel Size | Type of fitting | SD / Rivnut<br>Fitting distance | Weight | Item no.     |
|---------------------------------|------------|-----------------|---------------------------------|--------|--------------|
| Optimum H - Black hard anodized | 24         | 6 splices       | 510                             | 260    | 34244G (NEW) |
| Optimum H - Black hard anodized | 24         | 6 splices       | 515                             | 275    | 34244S       |
| Optimum H - Black hard anodized | 25         | 6 splices       | 530                             | 285    | 34254G       |
| Optimum H - Black hard anodized | 25         | 6 splices       | 535                             | 290    | 34254S       |
| Optimum H - Black hard anodized | 26         | 6 splices       | 565                             | 300    | 34264S       |
| Optimum H - Black hard anodized | 24         | 6 rivnuts       | 515                             | 270    | 32244R       |

