





# PU tyres for manual wheelchairs


## 12" – 36"




- » Puncture-proof and maintenance-free
- » No need to inflate
- » Constant tyre pressure
- » Suitable for indoor use with light manual wheelchairs
- » Moderately shock absorbing
- » Low rolling resistance

## Polyurethane tyres 12"–16"






Inches	ETRTO	 mm	 mm	 mm	Part no. Grey	Part no. Black
12½" x 2¼"	57-203	315	50	27	131262	131262B







Inches	ETRTO	 mm	 mm	 mm	Part no. Grey	Part no. Black
12½" x 2¼"	57-203	315	39	25,5	131253	131253B
12½" x 2¼"	57-203	315	42	27	131252	131252B
16" x 1¾"	47-305	400	44	27	131651	131651B
16" x 1¾"	47-305	400	44	27	131651H (hard)	131651H (hard)






Inches	ETRTO	 mm	 mm	 mm	Part no. Grey	Part no. Black
12½" x 2¼"	57-203	300	45	27	131251	131251B



Inches	ETRTO	 mm	 mm	 mm	Part no. Grey	Part no. Black
12½" x 2¼"	62-203		47	27	131261	131261B





Inches	ETRTO	 mm	 mm	 mm	Part no. Grey	Part no. Black
12½" x 2"	62-203	300	46	36	131200	131200B





Inches	ETRTO	 mm	 mm	Part no. Grey_	Part no. Black
12½" x 2¼"	62-203	56	36	131264	On request


## Polyurethane tyres 20"-36"





Tyre sizes		 mm	Part no. Grey	Part no. Black
Inches	ETRTO			
22" x 1"	25-489	19	132251	132251H
24" x 1"	25-540	19	132451	132451H
26" x 1"	25-590	19	132651	132651H




Tyre sizes		 mm	Part no. Grey	Part no. Black
Inches	ETRTO			
20" x 1 3/8"	37-451	19	132061	132061H
22" x 1 3/8"	37-489	19	132261	132261H
22" x 1 3/8"	37-498/501	19	132262	132262H
24" x 1 3/8"	37-540	19	132461	132461H
26" x 1 3/8"	37-590	19	132661	132661H




Tyre sizes			Part no. Grey	Part no. Black
Inches	ETRTO			
24" x 1 3/8"	37-540	19	132490	132465H





Tyre sizes			Part no. Grey	Part no. Black
Inches	ETRTO			
20" x 1"	25-451	19	132081	132081H
22" x 1"	25-489	19	132281	on request
24" x 1"	25-540	19	132481	132481H




Tyre sizes			Part no. Grey	Part no. Black
Inches	ETRTO			
26 x 1"	20-564	19	132600	on request
29" x 1"	26-683	19	132900	on request
32" x 1"	26-735	19	133200	on request
36" x 1"	26-839	19	133600	on request



Tyre sizes			Part no. Grey	Part no. Black
Inches	ETRTO			
20" x 1 3/8"	37-451	19	132071	132071H
22" x 1 3/8"	37-489	19	132271	132271H
22" x 1 3/8"	37-498/501	19	132272	132272H
24" x 1 3/8"	37-540	19	132471	132471H
26" x 1 3/8"	37-590	19	132671	on request



Tyre sizes			Part no. Grey	Part no. Black
Inches	ETRTO			
20" x 1 3/8"	37-451	19	132077	132077H
22" x 1 3/8"	37-489	19	132277	132277H
22" x 1 3/8"	37-489/501	19	132278	132277H
24" x 1"	23-540	19	132487	132487H
24" x 1 3/8"	37-540	19	132477	132477H