



Antitip system




Antitip – type A

- » Operated by foot
- » Can be fitted on most wheelchairs
- » Big adjustment range – stepless
- » Can swing away under the chair
- » Material: aluminum
- » Item reference no. 9311AP



Antitip – type B

- » Operated by hand or foot
- » Can be fitted inside a tube.
- » Can swing upwards when not used
- » Can be adjusted in several steps
- » Material: aluminum or steel
- » Item reference no. for aluminium version: left side MI010ZA, right side MI110ZA
- » Item reference no. for steel version: MI010ZA.









Antitip – type C

- » Hand operated
- » Can be fitted on tubes outer diam. Ø25mm
- » Big adjustment range
- » Material: Aluminum/Plastic
- » Item reference no.: 7212DE

Wheel for antitip systems



|  mm |  mm |  mm |  mm |  kg |  g | Part no. |
|--|--|--|--|--|---|----------|
| 47 | 10 | ø6 | 8/10 | - | - | 886015 |