

Back canes



Height adjustable back cane – TYPE 1

- » Stepless height adjustable range 120 mm
- » Release the knob and adjust the height by lifting up or pushing down
- » Available for different tube diameters and hole positions on request


Description	Reference part. no.
Height adjustable back cane	7212CK



Angle adjustable back cane – TYPE 2

- » Stepless angle adjustable range: 23 degrees
- » Starting angle is 85,5 degrees and end angle is 62,5 degrees
- » Release the knob and adjust the angle by lifting up or pushing down
- » Available for different tube diameters and hole positions on request

Description	Reference part. no.
Angle adjustable back cane	72111L



Angle adjustable back cane – TYPE 3

- » Angle adjustable range: 26 degrees
- » Starting angle is 90 degrees and can be adjusted to 77 degrees and 64 degrees
- » Press the button and adjust the angle by lifting up or pushing down
- » Available for different tube diameters and hole positions on request

Description	Reference part. no.
Height adjustable back cane	72111XA



Angle adjustable back cane – TYPE 4

- » Angle adjustable range: 40 degrees
- » Can be adjusted to following angle positions: 82,5 degrees; 72,5 degrees; 62,5 degrees; 52,5 degrees; 42,5 degrees
- » Pull out the position lock and adjust the angle by lifting up or pushing down
- » Available for different tube diameters and hole positions on request

Description	Reference part. no.
Angle adjustable back cane	72112X



Half foldable back cane – TYPE 5

- » Compact for transferring when half-folded
- » Produced in high strength casted steel
- » Press the handle and fold it. Lift up the handle to push back until you hear the "click" sound
- » Available for different tube diameters and hole positions on request

Description	Reference part. no.
Half foldable back cane	72112L