


Back hinges




Half foldable back cane hinge

- » Compact hinge for foldable back rests.
- » Produced in high strength casted steel.
- » Easy push the lever and turn upper part downward.
- » Fast change the back height for transportation and storage.
- » Available for different tube diameters and hole patterns on request.
- » Reference item no: 67315JM



Hinge for angle adjustable back rests

- » Newly designed and patented hinge for angle adjustable back rests.
- » Produced in high strength casted steel.
- » Incremental adjustment. 0°, -13°, -26°. Other angles available on request.
- » Easy and fast to change the angle of the back rest. Push the release button and adjust for the desired angle position. No tools needed.
- » Very rigid and stabile.
- » Available for different tube diameters and hole patterns on request.
- » Reference item no: LS0101



Hinge for half foldable back rests

- » Compact hinge for foldable back rests.
- » Produced in high strength casted steel.
- » Available for different tube diameters and hole patterns on request.
- » Reference item no.: 67315J