

# SIMI pushrims with anti-slip insert

The SIMI pushrim is an excellent choice for users with full or slightly reduced hand strength as it has a large surface and a shape that fits well into the hand. The height of the profile is 34 mm. On the top it has an anti-slip, yet not too sticky silicone insert made of environmentally friendly food grade silicone which provides an excellent grip.

- Large aluminium surface for easier braking
- The silicone insert is not too sticky which results in a reliable and comfortable grip as sand and other materials do not stick onto it.
- Good anti-slip performance even when wet.
- 6 splices (39210Q)
- Available in black hard anodized or shiny silver anodized, Silver is slightly more soft to touch but not hard anodized, so it scratches easier.
- Replaceable silicone insert



Wheel Size	SD fitting distance mm	Fit Omobic rim types	Weight g	Silver anodized (shiny) SIMI S	Black hard anodized SIMI H
24"	515	DW2 + DW4 + DW7	825	3424DB	3424DA
24"	523	SW1	825	3424DD	3424DC
25"	535	DW2 + DW4 + DW7	860	3425DB	3425DA
26"	565	DW2 + DW4 + DW7	900	3426DB	3426DA

