

HURRICANE 6 x 3

22"-25"

HURRICANE 6 x 3 is an optimal choice if you are looking for a modern and dynamic looking wheel. The 6 x 3 spokes combined with a 28 mm tall and narrow rim results in a lean and attractive look. HURRICANE is a strong yet light wheel and fitted with Omobic Speedlite tyre it will guarantee the user a very easy ride.





Black KTL
(DW5)



Aluminium pushrim
6 splices with two holes
Black or silver anodized
Ø19 mm



Black hub (F3) 53/64 mm
18 x Ø2,3mm stainless steel
spokes



Pushrim fitted with 6
round rivnuts



Max load
75 kg



Omobic SpeedLite
25 mm

Raw wheels

SIZE	Rivet Diameter mm	Weight g	Part no. Ø12 mm	Part no. Ø1/2"
22"/Ø489	469	740	232H31	232H81
24"/Ø540	520	790	234H31	234H81
25"/Ø559	550	830	235H31	235H81

Wheels with pushrim and tyre

SIZE	Pushrim	Weight g	Part no. Ø12 mm	Part no. Ø1/2"
24"/Ø540	Silver anodized aluminium	1530	264H31	264H81
24"/Ø540	Black anodized aluminium	1530	264H32	264H82

Quick-release axles in steel or titanium

Material	Weight (g) Ø12 mm / Ø1/2"	Length mm	Release button colour	Part. no. Ø12 mm	Part. no. Ø1/2"
Steel	110	108	Black	6180S1	6181S3
Titanium	75	108	Black	6180S9	6181S9