


BUDDY

24"

Our off-road wheel BUDDY will safely carry the user and prevent the wheelchair from sinking on soft grounds. The wide rim and fat tyre add style to any wheelchair. Our version with a 507 mm (22") wheel rim fitted with a 60 mm wide tyre has a diameter which is close to a standard 24" wheel and so the setup of the wheelchair remains the same.






Black
extra wide



Aluminium pushrim
6 splices with two holes
Black anodized
Ø19 mm



Black hub with high flanges
(F5) 53/64mm
36 stainless steel spokes,
2 mm



Big Apple 24"/507-
24" x 2" - other
DIRECTIONAL tyres
on request



Max load
60 kg

Raw wheels

SIZE	SD diameter mm	Weight g	Part no. Ø12 mm	Part no. Ø1/2"
24"/Ø507	480	1100	23F403	23F404

Wheels with pushrim and tyre

SIZE	Tyre DIRECTIONAL	Weight g	Part no. Ø12 mm	Part no. Ø1/2"
24"/Ø507	LEFT	2400	26F405L	26F406L
24"/Ø507	RIGHT	2400	26F405R	26F406R

Quick-release axles in steel or titanium

Material	Weight (g) Ø12 mm /Ø1/2"	Lenght mm	Release button colour	Part. no. Ø12 mm	Part. no. Ø1/2"
Steel	97/108	108	Black	618254	618354
Titanium	66/73	108	Black	61825A	61835A