

Disc-brake wheel with various options for brake handles 24"

Fast and easy release of the wheels. Simply push the button, lift up the brake system and release your wheel fitted with a QR axle! Available with 120 mm disc and equipped with a safety guard which prevents the rotor from touching the ground when the wheel is deattached. Chose between three options of brake handles for user or attendant activations. **Only available with 12 mm axle.**



Wheels with pushrim & tyre – without disc-brake



Please find item numbers for sets with brakes on the next pages!

SIZE	SD diameter mm	Hub type	No. of spokes	 g	Part no. ø12mm
24"/Ø540	523	69650V	18	2330	234AAZ



Black rim (SW1)



Silver anodized aluminium pushrim with 6 splices and two holes. Ø19 mm



Black hub 70/48 mm
18 x ø2,6 mm black KTL spokes



Pushrim fitting with 6 splices, screw and lock-nut.



Max load for the system
150 kg



Omobic SpeedLite 25 mm

Complete wheels with brake handle – engaged by the attendant



- Disc-brake system with Ø120 mm rotor
- 12 or 12,77 mm QR axle
- Mounting clamps included (left or right side)
- NB! Item no is for ONE wheel
- Brake handle 69DHD1 (left) 69DHD2 (right) included

SIZE	AXLE	BRAKE HANDLE	Weight g	Part no. ONE wheel RIGHT side	Part no. ONE wheel LEFT side
24"/Ø540	12 mm QR axle	69DHD1 (left)	2330	69D02R	69D02L
24"/Ø540	12,7 mm QR axle	69DHD2 (right)	2330	69D06R	69D06L

Complete wheels with one cable in each brake handle – engaged by the user



- Disc-brake system with Ø120 mm rotor
- 12 mm QR axle
- Including brake handle 67103M
- Mounting clamps included (left or right side)
- NB! Item no is for ONE wheel
- Brake handle 69DHD1 (left) 69DHD2 (right) included


SIZE	AXLE	Weight g	Part no. ONE wheel RIGHT side	Part no. ONE wheel LEFT side
24"/Ø540	12 mm QR axle	2330	69D08R	69D08L

Set of wheels with two cables in 1 brake handle – engaged by the user



- Disc-brake system with Ø120 mm rotor
- 12 mm QR axle
- Including brake handles 67103M (left side) and 67103N (right side)
- Mounting clamps included (left or right side)
- NB! Item no is for TWO wheel

SIZE	AXLE	BRAKE HANDLE	Weight g	Part no. TWO wheels Mounting clamp for han- dle RIGHT side	Part no. TWO wheels Mounting clamp for han- dle LEFT side
24"/Ø540	12 mm QR axle	67103M (left) and 67103N (right)	2330	69D10R	69D10L




Position 1: Parking brake engaged
Position 2: Parking brake not engaged



2 cable outlets



The brake handle is assembled to the chair frame by using clamps as shown above.



The wheel is detachable by pressing the quick-release axle.



The safety guard prevents the rotor from touching the ground.