

Product Technical Data Sheet



Product application

Lock A3 wheel lock is dedicated for active wheelchair users to be used for stabilization and safety when the wheelchair is stationary. Typical usage: transfers in and out of the wheelchair, parking, sports and active use, safety on slopes on uneven ground, transport and vehicle loading.

Key features and USPs

- Easy operation – can be used also by users with reduced hand strength.
- Thanks to an optimized design, reduces the risk of hitting users' fingers.
- Collects less dirt from the tires (more "closed" when unlocked)
- Fits most chairs and it is suitable for all types of wheelchairs

Product specifications

- material: aluminum
- surface treatment: black anodized
- weight: 115 g
- 90 degree
- different assemble hole options

Options and item numbers

Description	Weight (g)	Part number
Active wheel lock with right hand	115	67303Y-DA
Active wheel lock with left hand	115	67303Z-DA

Mounting clamps for different frame dimensions

Fitting tube dimension	Part number
Ø19/Ø20mm	674005
Ø19/Ø22mm	674004
Ø19/Ø25mm	674001
Ø19/Ø30mm	674003