

Airy handrim with anti-slip insert

Grip well. Weight less.

The AIRY handrim is the lightest ergonomic handrim available on the market with good grip you were looking for. The profile is optimized for users with small and big hands and still offering very ergonomic design.

- Good grip even when wet – thanks to the anti-slip silicone insert and optimal surface.
- Optimized profile design – designed large area without silicone strip – so you don't burn your hands during braking.
- Easy to keep clean – because the insert is made of silicone that is not too sticky and doesn't catch sand or dirt from the ground.
- Low weight – comparing to other handrims with silicone strip on the market.
- Fits well into hand – ergonomic shape.



Wheel Size	SD fitting distance mm	Colour	Weight g	Surface	Item no
24"	515	Black	560	Polished	3424HS
24"	515	Silver	560	Polished	3424HS-S
24"	515	Black	560	Brushed	3424HB
24"	515	Silver	560	Brushed	3424HB-S
24"	510	Black	560	Polished	3424HA
24"	510	Silver	560	Polished	3424HA-S
24"	510	Black	560	Brushed	3424HC
24"	510	Silver	560	Brushed	3424HC-S