

Strenght. Spoke by spoke.

The REBEL is extremely strong wheel with fibre spokes for demanding whelchair users that can be used for everyday life and with power add-ons. It is a wheel for active users who are looking for a good -looking wheel with a very good trade -off between weight, performance and durability.

- Tested for 180 kg also with cambering.
- Fibre spokes - 50% ligther than steel spokes and 3 times stronger tahn steel spokes. Absorbent. Durable. Maneuverable.
- Reliable structure, suitable for wheelchair users with higher body weight .
- Excellent shock-absorption eases road-surface vibrations, engineered for users pursuing a superior riding experience.
- Highly durable spokes.
- Affortable price.



Wheel Size	SD fitting distance mm	Colour	Weight g	Part no. ø12 mm	Part no. ø1/2"
24"	510	Black	840	F240000002	F240000001
25"	530	Black	930	F250000002	F250000001

Wheel Size	Spoke item no
24"	22B1F243CN
25"	22B1F252